

**pHBuffer Solution Recal mixed pack of 2xpH4 7 & 10 ±0.01 @20°C**

Art. ID	REA-MXC60
Unit	6 x 90 mL
Deliverydetails	No Dangerous Good /not restricted

**Description**

Supplied in convenient and economical wide mouth, single-use, 6x90ml container designed for Field and Lab direct calibration. Tested at 20 or 25°C and Certified in Reagecon ISO 17025 Accredited Laboratory. Traceable to a NIST Standard. Accuracy ±0.04