

**pHBuffer Solution pH4.00 (Red)  $\pm 0.01$  @25°C (RECAL)**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | REA-04C65                         |
| Unit            | 6 x 90 mL                         |
| Deliverydetails | No Dangerous Good /not restricted |

**Description**

Supplied in convenient and economical wide mouth, single-use, 6x90ml container designed for Field and Lab direct calibration. Tested at 20 or 25°C and Certified in Reagecon ISO 17025 Accredited Laboratory. Traceable to a NIST Standard. Accuracy  $\pm 0.05$