

### **Conductivity Standard 1413 Microsiemens/cm @25°C**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | REA-CSKCL                         |
| Unit            | 500 mL                            |
| Deliverydetails | No Dangerous Good /not restricted |

#### **Description**

Conductivity Standard Values @25°C, Supplied with ISO17025 Certificate and Directly Traceable to NIST. Actual value, expiry and temperature dependance data on label. All Reagecon Conductivity standards are aqueous based and non-hazardous.