

CEN PCB Congener Mix-1

Art. ID SAF-47927
Unit 1 mL
Deliverydetails De Minimis UN1208 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',5-Trichlorobipheny I (PCB 18)	[37680-65-2]	10	µg/mL		
	2,4,4'-Trichlorobipheny I (PCB 28)	[7012-37-5]	10	µg/mL		
	2,4',5-Trichlorobipheny I (PCB 31)	[16606-02-3]	10	µg/mL		
	2,2',3,5'-Tetrachlorobi phenyl (PCB 44)	[41464-39-5]	10	µg/mL		
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	10	µg/mL		
	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	10	µg/mL		
	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	10	µg/mL		
	2,2',3,4,4',5'-Hexachlo robiphenyl (PCB 138)	[35065-28-2]	10	µg/mL		
	2,2',3,4',5',6-Hexachlo robiphenyl (PCB 149)	[38380-04-0]	10	µg/mL		
	2,2',4,4',5,5'-Hexachlo robiphenyl (PCB 153)	[35065-27-1]	10	µg/mL		
	2,2',3,4,4',5,5'-Heptac hlorobiphenyl (PCB 180)	[35065-29-3]	10	µg/mL		
	2,2',3,3',4,4',5,5'-Oct achlorobiphenyl (PCB 19 4)	[35694-08-7]	10	µg/mL		