

**Tuning Solution for Agilent ICP-MS 7500cs, contains 1 µg/L each of Li, Mg, Y, Ce, Tl and Co in a matrix of 2% HNO<sub>3</sub>**

Art. ID LM24-70017-ICP-MS-250mL  
Unit 250 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	1	µg/L		
	Magnesium (Mg)	[7439-95-4]	1	µg/L		
	Yttrium (Y)	[7440-65-5]	1	µg/L		
	Cerium (Ce)	[7440-45-1]	1	µg/L		
	Thallium (Tl)	[7440-28-0]	1	µg/L		
	Cobalt (Co)	[7440-48-4]	1	µg/L		