

## **Conductivity Standard 20,000 Microsiemens/cm @25°C**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | REA-CSKC20M                       |
| Unit            | 500 mL                            |
| Deliverydetails | No Dangerous Good /not restricted |

### Description

Conductivity Premium Values @25°C, Supplied with ISO17025 Certificate and Directly Traceable to NIST. Actual value, expiry and temperature dependance data on label. All Reagecon Conductivity standards are aqueous based and non-hazardous.