

Multi-Element Standard, 5 components each 1000 mg/L in 5% HNO₃

Art. ID LM24-SL-70039
Unit 100 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	1000	mg/L		
	Germanium (Ge)	[7440-56-4]	1000	mg/L		
	Lithium (Li)	[7439-93-2]	1000	mg/L		
	Scandium (Sc)	[7440-20-2]	1000	mg/L		
	Indium (In)	[7440-74-6]	1000	mg/L		