

ICP Multi element mixture 20 components in 5% HNO₃

Art. ID LM24-SL-70041
Unit 100 ml
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Arsenic (As) | [7440-38-2] | 100 | µg/mL | | |
| | Boron (B) | [7440-42-8] | 100 | µg/mL | | |
| | Cadmium (Cd) | [7440-43-9] | 100 | µg/mL | | |
| | Calcium (Ca) | [7440-70-2] | 100 | µg/mL | | |
| | Copper (Cu) | [7440-50-8] | 100 | µg/mL | | |
| | Manganese (Mn) | [7439-96-5] | 100 | µg/mL | | |
| | Selenium (Se) | [7782-49-2] | 100 | µg/mL | | |
| | Silver (Ag) | [7440-22-4] | 100 | µg/mL | | |
| | Potassium (K) | [7440-09-7] | 100 | µg/mL | | |
| | Sodium (Na) | [7440-23-5] | 100 | µg/mL | | |
| | Phosphorus (P) | [7723-14-0] | 100 | µg/mL | | |
| | Aluminium (Al) | [7429-90-5] | 100 | µg/mL | | |
| | Chromium (Cr) | [7440-47-3] | 100 | µg/mL | | |
| | Tin (Sn) | [7440-31-5] | 100 | µg/mL | | |
| | Zinc (Zn) | [7440-66-6] | 100 | µg/mL | | |
| | Iron (Fe) | [7439-89-6] | 100 | µg/mL | | |
| | Lead (Pb) | [7439-92-1] | 100 | µg/mL | | |
| | Magnesium (Mg) | [7439-95-4] | 100 | µg/mL | | |
| | Nickel (Ni) | [7440-02-0] | 100 | µg/mL | | |
| | Mercury (Hg) | [7439-97-6] | 100 | µg/mL | | |