

### Custom mixture, 13 components in Methanol

Art. ID LM24-SL-70044  
Unit 5 x 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,1,1-Trichloroethane	[71-55-6]	200	µg/mL		
	1,1,2,2-Tetrachloroethane	[79-34-5]	200	µg/mL		
	1,1,2-Trichloroethane	[79-00-5]	200	µg/mL		
	1,1-Dichloroethene	[75-35-4]	200	µg/mL		
	1,2-Dichloroethane	[107-06-2]	200	µg/mL		
	1,2-Dichloropropane	[78-87-5]	200	µg/mL		
	Carbon tetrachloride (Tetrachloromethane)	[56-23-5]	200	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	200	µg/mL		
	cis-1,2-Dichloroethene	[156-59-2]	200	µg/mL		
	Tetrachloroethene	[127-18-4]	200	µg/mL		
	trans-1,2-Dichloroethene	[156-60-5]	200	µg/mL		
	Trichloroethene	[79-01-6]	200	µg/mL		
	Vinylchloride	[75-01-4]	200	µg/mL		