

Pinnatoxin-G (PnTX-G) in Methanol

Art. ID	NRC-CRM-PnTX-G
Unit	0,5 mL
Deliverydetails	De Minimis UN1230 (6.1) II E2

Description

NRC-CRM-PnTX-G is a certified calibration solution for the determination of pinnatoxin-G (PnTX-G) by chemical analysis methods. Each ampoule contains ~0.5 mL of a 2.77 µmol/L solution of PnTX-G in methanol (with 0.01% acetic acid).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Pinnatoxin-G (PnTX-G)		2,77	µmol/L		