

**Mixture Ba + Si, aqueous calibration solution in 2% HNO<sub>3</sub>**

Art. ID LM24-CUS-70052  
Unit 250 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	5	mg/L		
	Silicon (Si)	[7440-21-3]	20	mg/L		