

PAH Mix 4, Varied conc, in Acetonitrile

Art. ID ACCU-AE-00033
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzo(b)fluoranthene	[205-99-2]	20	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	20	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	20	µg/mL		
	Benzo(a)pyrene	[50-32-8]	20	µg/mL		
	Fluoranthene	[206-44-0]	50	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	40	µg/mL		
	Perylene	[198-55-0]	20	µg/mL		