

Multi-element Calibration Standard #3, 10 µg/mL in Water, % Nitric acid, % Hydrofluoric acid

Art. ID ACCU-AG-MECAL4-ASL-R1-5
Unit 500 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	10	µg/mL		
	Germanium (Ge)	[7440-56-4]	10	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	10	µg/mL		
	Niobium (Nb)	[7440-03-1]	10	µg/mL		
	Phosphorus (P)	[7723-14-0]	10	µg/mL		
	Rhenium (Re)	[7440-15-5]	10	µg/mL		
	Sulfur	[7704-34-9]	10	µg/mL		
	Silicon (Si)	[7440-21-3]	10	µg/mL		
	Tin (Sn)	[7440-31-5]	10	µg/mL		
	Tantalum (Ta)	[7440-25-7]	10	µg/mL		
	Titanium (Ti)	[7440-32-6]	10	µg/mL		
	Tungsten (W)	[7440-33-7]	10	µg/mL		
	Zirconium (Zr)	[7440-67-7]	10	µg/mL		