

7500 Series PA Tuning 1, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-AG-TUN1-ASL-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	20	µg/mL		
	Beryllium (Be)	[7440-41-7]	20	µg/mL		
	Cadmium (Cd)	[7440-43-9]	20	µg/mL		
	Arsenic (As)	[7440-38-2]	20	µg/mL		
	Nickel (Ni)	[7440-02-0]	10	µg/mL		
	Lead (Pb)	[7439-92-1]	10	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	5	µg/mL		
	Sodium (Na)	[7440-23-5]	5	µg/mL		
	Aluminium (Al)	[7429-90-5]	5	µg/mL		
	Uranium (U)	[7440-61-1]	5	µg/mL		
	Copper (Cu)	[7440-50-8]	5	µg/mL		
	Thorium (Th)	[7440-29-1]	5	µg/mL		
	Barium (Ba)	[7440-39-3]	5	µg/mL		
	Cobalt (Co)	[7440-48-4]	5	µg/mL		
	Strontium (Sr)	[7440-24-6]	5	µg/mL		
	Vanadium (V)	[7440-62-2]	5	µg/mL		
	Chromium (Cr)	[7440-47-3]	5	µg/mL		
	Manganese (Mn)	[7439-96-5]	5	µg/mL		
	Lithium-6 (LI-6)	[14258-72-1]	5	µg/mL		
	Scandium (Sc)	[7440-20-2]	5	µg/mL		
	Indium (In)	[7440-74-6]	5	µg/mL		
	Lutetium (Lu)	[7439-94-3]	5	µg/mL		
	Bismuth (Bi)	[7440-69-9]	5	µg/mL		
	Yttrium (Y)	[7440-65-5]	2,5	µg/mL		
	Ytterbium (Yb)	[7440-64-4]	2,5	µg/mL		