

**7500 Series PA Tuning 2, Varied conc, in 10% Hydrochloric acid, 1% Nitric acid, tr Hydrofluoric acid**

Art. ID ACCU-AG-TUN2-ASL-5  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	10	µg/mL		
	Antimony (Sb)	[7440-36-0]	10	µg/mL		
	Tin (Sn)	[7440-31-5]	10	µg/mL		
	Germanium (Ge)	[7440-56-4]	10	µg/mL		
	Ruthenium (Ru)	[7440-18-8]	10	µg/mL		
	Palladium (Pd)	[7440-05-3]	10	µg/mL		
	Titanium (Ti)	[7440-32-6]	5	µg/mL		
	Iridium (Ir)	[7439-88-5]	5	µg/mL		