

## **7500 Series PA Tuning Solution Set**

Art. ID	ACCU-AG-TUN-ASL-1-SET
Unit	2 x 100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-AG-TUN1-ASL-1	Aluminium (Al)	[7429-90-5]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Arsenic (As)	[7440-38-2]	20	µg/mL		
ACCU-AG-TUN1-ASL-1	Barium (Ba)	[7440-39-3]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Beryllium (Be)	[7440-41-7]	20	µg/mL		
ACCU-AG-TUN1-ASL-1	Bismuth (Bi)	[7440-69-9]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Cadmium (Cd)	[7440-43-9]	20	µg/mL		
ACCU-AG-TUN1-ASL-1	Chromium (Cr)	[7440-47-3]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Cobalt (Co)	[7440-48-4]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Copper (Cu)	[7440-50-8]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Indium (In)	[7440-74-6]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Lead (Pb)	[7439-92-1]	10	µg/mL		
ACCU-AG-TUN1-ASL-1	Lithium-6 (LI-6)	[14258-72-1]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Lutetium (Lu)	[7439-94-3]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
ACCU-AG-TUN1-ASL-1	Manganese (Mn)	[7439-96-5]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Nickel (Ni)	[7440-02-0]	10	µg/mL		
ACCU-AG-TUN1-ASL-1	Scandium (Sc)	[7440-20-2]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Sodium (Na)	[7440-23-5]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Strontium (Sr)	[7440-24-6]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Thallium (Tl)	[7440-28-0]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Thorium (Th)	[7440-29-1]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Uranium (U)	[7440-61-1]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Vanadium (V)	[7440-62-2]	5	µg/mL		
ACCU-AG-TUN1-ASL-1	Ytterbium (Yb)	[7440-64-4]	2,5	µg/mL		
ACCU-AG-TUN1-ASL-1	Yttrium (Y)	[7440-65-5]	2,5	µg/mL		
ACCU-AG-TUN1-ASL-1	Zinc (Zn)	[7440-66-6]	20	µg/mL		
ACCU-AG-TUN2-ASL-1	Antimony (Sb)	[7440-36-0]	10	µg/mL		
ACCU-AG-TUN2-ASL-1	Germanium (Ge)	[7440-56-4]	10	µg/mL		
ACCU-AG-TUN2-ASL-1	Iridium (Ir)	[7439-88-5]	5	µg/mL		
ACCU-AG-TUN2-ASL-1	Molybdenum (Mo)	[7439-98-7]	10	µg/mL		
ACCU-AG-TUN2-ASL-1	Palladium (Pd)	[7440-05-3]	10	µg/mL		

ACCU-AG-TUN2-ASL-1	Ruthenium (Ru)	[7440-18-8]	10	µg/mL
ACCU-AG-TUN2-ASL-1	Tin (Sn)	[7440-31-5]	10	µg/mL
ACCU-AG-TUN2-ASL-1	Titanium (Ti)	[7440-32-6]	5	µg/mL