

## **Sulfur & Metals in Residual Fuel Oil**

Art. ID                    ACCU-ASTM-P-0103-SET  
Unit                      12 x 100 mL  
Deliverydetails        No Dangerous Good

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-ASTM-P-0103-01	Iron (Fe)	[7439-89-6]	Varied			
ACCU-ASTM-P-0103-01	Nickel (Ni)	[7440-02-0]	Varied			
ACCU-ASTM-P-0103-01	Sulfur (S)	[7704-34-9]	Varied			
ACCU-ASTM-P-0103-01	Vanadium (V)	[7440-62-2]	Varied			
ACCU-ASTM-P-0103-02	Iron (Fe)	[7439-89-6]	300 Varied			
ACCU-ASTM-P-0103-02	Nickel (Ni)	[7440-02-0]	10 Varied			
ACCU-ASTM-P-0103-02	Sulfur (S)	[7704-34-9]	0,5 Varied			
ACCU-ASTM-P-0103-02	Vanadium (V)	[7440-62-2]	500 Varied			
ACCU-ASTM-P-0103-03	Iron (Fe)	[7439-89-6]	500 Varied			
ACCU-ASTM-P-0103-03	Nickel (Ni)	[7440-02-0]	100 Varied			
ACCU-ASTM-P-0103-03	Sulfur (S)	[7704-34-9]	1 Varied			
ACCU-ASTM-P-0103-03	Vanadium (V)	[7440-62-2]	25 Varied			
ACCU-ASTM-P-0103-04	Iron (Fe)	[7439-89-6]	100 Varied			
ACCU-ASTM-P-0103-04	Nickel (Ni)	[7440-02-0]	80 Varied			
ACCU-ASTM-P-0103-04	Sulfur (S)	[7704-34-9]	Varied			
ACCU-ASTM-P-0103-04	Vanadium (V)	[7440-62-2]	250 Varied			
ACCU-ASTM-P-0103-05	Iron (Fe)	[7439-89-6]	200 Varied			
ACCU-ASTM-P-0103-05	Nickel (Ni)	[7440-02-0]	40 Varied			
ACCU-ASTM-P-0103-05	Sulfur (S)	[7704-34-9]	2 Varied			
ACCU-ASTM-P-0103-05	Vanadium (V)	[7440-62-2]	100 Varied			
ACCU-ASTM-P-0103-06	Iron (Fe)	[7439-89-6]	400 Varied			
ACCU-ASTM-P-0103-06	Nickel (Ni)	[7440-02-0]	5 Varied			
ACCU-ASTM-P-0103-06	Sulfur (S)	[7704-34-9]	2,5 Varied			
ACCU-ASTM-P-0103-06	Vanadium (V)	[7440-62-2]	400 Varied			
ACCU-ASTM-P-0103-07	Iron (Fe)	[7439-89-6]	Varied			
ACCU-ASTM-P-0103-07	Nickel (Ni)	[7440-02-0]	60 Varied			
ACCU-ASTM-P-0103-07	Sulfur (S)	[7704-34-9]	3 Varied			
ACCU-ASTM-P-0103-07	Vanadium (V)	[7440-62-2]	300 Varied			
ACCU-ASTM-P-0103-08	Iron (Fe)	[7439-89-6]	500 Varied			
ACCU-ASTM-P-0103-08	Nickel (Ni)	[7440-02-0]	Varied			
ACCU-ASTM-P-0103-08	Sulfur (S)	[7704-34-9]	3,5 Varied			

ACCU-ASTM-P-0103-08	Vanadium (V)	[7440-62-2]	200 Varied
ACCU-ASTM-P-0103-09	Iron (Fe)	[7439-89-6]	100 Varied
ACCU-ASTM-P-0103-09	Nickel (Ni)	[7440-02-0]	100 Varied
ACCU-ASTM-P-0103-09	Sulfur (S)	[7704-34-9]	Varied
ACCU-ASTM-P-0103-09	Vanadium (V)	[7440-62-2]	Varied
ACCU-ASTM-P-0103-10	Iron (Fe)	[7439-89-6]	300 Varied
ACCU-ASTM-P-0103-10	Nickel (Ni)	[7440-02-0]	50 Varied
ACCU-ASTM-P-0103-10	Sulfur (S)	[7704-34-9]	4,5 Varied
ACCU-ASTM-P-0103-10	Vanadium (V)	[7440-62-2]	250 Varied
ACCU-ASTM-P-0103-11	Iron (Fe)	[7439-89-6]	200 Varied
ACCU-ASTM-P-0103-11	Nickel (Ni)	[7440-02-0]	20 Varied
ACCU-ASTM-P-0103-11	Sulfur (S)	[7704-34-9]	5 Varied
ACCU-ASTM-P-0103-11	Vanadium (V)	[7440-62-2]	500 Varied
ACCU-ASTM-P-0103-12	Iron (Fe)	[7439-89-6]	50 Varied
ACCU-ASTM-P-0103-12	Nickel (Ni)	[7440-02-0]	100 Varied
ACCU-ASTM-P-0103-12	Sulfur (S)	[7704-34-9]	5,5 Varied
ACCU-ASTM-P-0103-12	Vanadium (V)	[7440-62-2]	50 Varied