

Vanadium Standard - High Range Set

| | |
|-----------------|---|
| Art. ID | ACCU-ASTM-P-0105-SET |
| Unit | 7 x 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1 |

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|---------------------|-------------------|-------------|---------------------|------|--------|--------|
| ACCU-ASTM-P-0105-01 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-01 | Vanadium (V) | [7440-62-2] | 0 | %w/w | | |
| ACCU-ASTM-P-0105-02 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-02 | Vanadium (V) | [7440-62-2] | 0,03 | %w/w | | |
| ACCU-ASTM-P-0105-03 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-03 | Vanadium (V) | [7440-62-2] | 0,04 | %w/w | | |
| ACCU-ASTM-P-0105-04 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-04 | Vanadium (V) | [7440-62-2] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-05 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-05 | Vanadium (V) | [7440-62-2] | 0,06 | %w/w | | |
| ACCU-ASTM-P-0105-06 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-06 | Vanadium (V) | [7440-62-2] | 0,08 | %w/w | | |
| ACCU-ASTM-P-0105-07 | Manganese (Mn) | [7439-96-5] | 0,05 | %w/w | | |
| ACCU-ASTM-P-0105-07 | Vanadium (V) | [7440-62-2] | 0,1 | %w/w | | |