

4 Component Deuterated ISTD Standard, Varied mass% in Isooctane

Art. ID ACCU-ASTM-P-0140-IS2
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene D6	[1076-43-3]	2	g/100g		
	Ethylbenzene-D10	[25837-05-2]	2	g/100g		
	Naphthalene D8	[1146-65-2]	1	g/100g		
	Toluene D8	[2037-26-5]	7	g/100g		