

### **CLP Surrogate Spiking Solution, Varied conc, in Methanol**

Art. ID ACCU-CLP-031-R-WL-25ML  
Unit 25 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Chlorophenol-3,4,5,6-D4	[93951-73-6]	150	µg/mL		
	1,2-Dichlorobenzene D4	[2199-69-1]	100	µg/mL		
	2-Fluorobiphenyl	[321-60-8]	100	µg/mL		
	2-Fluorophenol	[367-12-4]	150	µg/mL		
	Nitrobenzene D5	[4165-60-0]	100	µg/mL		
	Phenol D5	[4165-62-2]	150	µg/mL		
	p-Terphenyl D14	[1718-51-0]	100	µg/mL		
	2,4,6-Tribromophenol	[118-79-6]	150	µg/mL		