

Palladium Matrix Modifier Pd - 10.0 +/- 0.2 g/L in diluted HNO3 (approx. 15 %)

Art. ID LM24-CUS-70275
Unit 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Palladium (Pd) | [7440-05-3] | 10 ± 0,2 | g/L | | |