

### **CLP Calibration Std #3, Varied conc, in 2-5% Nitric Acid**

Art. ID ACCU-CLP-CAL-03-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	2000	µg/mL		
	Barium (Ba)	[7440-39-3]	2000	µg/mL		
	Beryllium (Be)	[7440-41-7]	50	µg/mL		
	Cobalt (Co)	[7440-48-4]	500	µg/mL		
	Copper (Cu)	[7440-50-8]	250	µg/mL		
	Iron (Fe)	[7439-89-6]	1000	µg/mL		
	Vanadium (V)	[7440-62-2]	500	µg/mL		