

CLP Initial Calibration Verification Std, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-CLP-ICV-01-R-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 200 | µg/mL | | |
| | Antimony (Sb) | [7440-36-0] | 60 | µg/mL | | |
| | Arsenic (As) | [7440-38-2] | 10 | µg/mL | | |
| | Barium (Ba) | [7440-39-3] | 200 | µg/mL | | |
| | Beryllium (Be) | [7440-41-7] | 5 | µg/mL | | |
| | Cadmium (Cd) | [7440-43-9] | 5 | µg/mL | | |
| | Calcium (Ca) | [7440-70-2] | 500 | µg/mL | | |
| | Chromium (Cr) | [7440-47-3] | 10 | µg/mL | | |
| | Cobalt (Co) | [7440-48-4] | 50 | µg/mL | | |
| | Copper (Cu) | [7440-50-8] | 25 | µg/mL | | |
| | Iron (Fe) | [7439-89-6] | 100 | µg/mL | | |
| | Lead (Pb) | [7439-92-1] | 5 | µg/mL | | |
| | Magnesium (Mg) | [7439-95-4] | 500 | µg/mL | | |
| | Manganese (Mn) | [7439-96-5] | 15 | µg/mL | | |
| | Nickel (Ni) | [7440-02-0] | 40 | µg/mL | | |
| | Potassium (K) | [7440-09-7] | 500 | µg/mL | | |
| | Selenium (Se) | [7782-49-2] | 5 | µg/mL | | |
| | Silver (Ag) | [7440-22-4] | 10 | µg/mL | | |
| | Sodium (Na) | [7440-23-5] | 500 | µg/mL | | |
| | Thallium (Tl) | [7440-28-0] | 10 | µg/mL | | |
| | Vanadium (V) | [7440-62-2] | 50 | µg/mL | | |
| | Zinc (Zn) | [7440-66-6] | 20 | µg/mL | | |