

Laboratory Control Sample - Internal Standard Mix, 2500 µg/mL in Methanol

Art. ID ACCU-CLP-LC-IS-100X-PAK
Unit 5 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,4-Dichlorobenzene D4	[3855-82-1]	2500	µg/mL		
	Chlorobenzene D5	[3114-55-4]	2500	µg/mL		
	1,4-Difluorobenzene	[540-36-3]	2500	µg/mL		