

**CLP Alternate Analytes Standard, Varied conc, in 2-5% Nitric Acid, tr Hydrofluoric acid**

Art. ID ACCU-CLP-PAN-02-5  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 100                 | µg/mL |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 100                 | µg/mL |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 100                 | µg/mL |        |        |
|                  | Boron (B)         | [7440-42-8] | 100                 | µg/mL |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 10                  | µg/mL |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 10                  | µg/mL |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 10                  | µg/mL |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 100                 | µg/mL |        |        |
|                  | Selenium (Se)     | [7782-49-2] | 100                 | µg/mL |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 10                  | µg/mL |        |        |
|                  | Sodium (Na)       | [7440-23-5] | 100                 | µg/mL |        |        |
|                  | Thallium (Tl)     | [7440-28-0] | 100                 | µg/mL |        |        |