

PCB Congeners Mix # 4 (NON Inducer Type), 10 µg/mL in Isooctane

Art. ID ACCU-C-SCA-04
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',5-Trichlorobipheny I (PCB 18)	[37680-65-2]	10	µg/mL		
	2,2',3,5'-Tetrachlorobi phenyl (PCB 44)	[41464-39-5]	10	µg/mL		
	2,2',4,5'-Tetrachlorobi phenyl (PCB 49)	[41464-40-8]	10	µg/mL		
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	10	µg/mL		
	2,3',4',5-Tetrachlorobi phenyl (PCB 70)	[32598-11-1]	10	µg/mL		
	2,4,4',5-Tetrachlorobip henyl (PCB 74)	[32690-93-0]	10	µg/mL		
	2,2',3,5,5',6-Hexachlor obiphenyl (PCB 151)	[52663-63-5]	10	µg/mL		
	2,2',3,3',4',5,6-Heptac hlorobiphenyl (PCB 177)	[52663-70-4]	10	µg/mL		
	2,2',3,4',5,5',6-Heptac hlorobiphenyl (PCB 187)	[52663-68-0]	10	µg/mL		
	2,2',3,3',4',5,5',6-Oct achlorobiphenyl (IUPAC PCB 199, PCB 201)	[52663-75-9]	10	µg/mL		