

## **Aromatics in Gasoline by GC/TC**

Art. ID ACCU-D-4420-CAL-SET  
Unit 7 x 1 mL  
Deliverydetails De Minimis

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-D-4420-CAL-01	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	64,45	%v/v		
ACCU-D-4420-CAL-01	Benzene	[71-43-2]	0,05	%v/v		
ACCU-D-4420-CAL-01	n-Butylbenzene	[104-51-8]	30	%v/v		
ACCU-D-4420-CAL-01	Toluene	[108-88-3]	0,5	%v/v		
ACCU-D-4420-CAL-01	Xylenes	[1330-20-7]	5	%v/v		
ACCU-D-4420-CAL-02	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	63,9	%v/v		
ACCU-D-4420-CAL-02	Benzene	[71-43-2]	0,1	%v/v		
ACCU-D-4420-CAL-02	n-Butylbenzene	[104-51-8]	25	%v/v		
ACCU-D-4420-CAL-02	Toluene	[108-88-3]	1	%v/v		
ACCU-D-4420-CAL-02	Xylenes	[1330-20-7]	10	%v/v		
ACCU-D-4420-CAL-03	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	62,25	%v/v		
ACCU-D-4420-CAL-03	Benzene	[71-43-2]	0,25	%v/v		
ACCU-D-4420-CAL-03	n-Butylbenzene	[104-51-8]	20	%v/v		
ACCU-D-4420-CAL-03	Toluene	[108-88-3]	2,5	%v/v		
ACCU-D-4420-CAL-03	Xylenes	[1330-20-7]	15	%v/v		
ACCU-D-4420-CAL-04	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	64,25	%v/v		
ACCU-D-4420-CAL-04	Benzene	[71-43-2]	0,75	%v/v		
ACCU-D-4420-CAL-04	n-Butylbenzene	[104-51-8]	10	%v/v		
ACCU-D-4420-CAL-04	Toluene	[108-88-3]	5	%v/v		
ACCU-D-4420-CAL-04	Xylenes	[1330-20-7]	20	%v/v		
ACCU-D-4420-CAL-05	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	58,75	%v/v		
ACCU-D-4420-CAL-05	Benzene	[71-43-2]	1,25	%v/v		
ACCU-D-4420-CAL-05	n-Butylbenzene	[104-51-8]	5	%v/v		
ACCU-D-4420-CAL-05	Toluene	[108-88-3]	10	%v/v		
ACCU-D-4420-CAL-05	Xylenes	[1330-20-7]	25	%v/v		
ACCU-D-4420-CAL-06	2,2,4-Trimethylpentane	[540-84-1]	66,5	%v/v		

(Isooctane)				
ACCU-D-4420-CAL-06	Benzene	[71-43-2]	2,5	%v/v
ACCU-D-4420-CAL-06	n-Butylbenzene	[104-51-8]	15	%v/v
ACCU-D-4420-CAL-06	Toluene	[108-88-3]	15	%v/v
ACCU-D-4420-CAL-06	Xylenes	[1330-20-7]	1	%v/v
ACCU-D-4420-CAL-07	2,2,4-Trimethylpentane	[540-84-1]	64,5	%v/v
(Isooctane)				
ACCU-D-4420-CAL-07	Benzene	[71-43-2]	5	%v/v
ACCU-D-4420-CAL-07	n-Butylbenzene	[104-51-8]	2,5	%v/v
ACCU-D-4420-CAL-07	Toluene	[108-88-3]	25	%v/v
ACCU-D-4420-CAL-07	Xylenes	[1330-20-7]	3	%v/v