

O-Acetylsalicyloyl Chloride

Art. ID TRC-A187935-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3261 II E2

Description

Application Note: O-Acetylsalicyloyl chloride is used in the derivatisation of amines, and also improve electrochemical detection of these amines when they are run through high-performance liquid chromatography. O-Acetylsalicyloyl chloride is also used as a reagent to synthesize benzimidazole derivatives, compounds that inhibit leukotriene synthesis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	O-Acetylsalicyloyl chloride	[5538-51-2]				