

ASTM/EPA Gasoline Refinery Standards, Varied conc, by weight

Art. ID ACCU-D-5769-QC-10ML
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Hexane	[110-54-3]	12	%		
	n-Heptane	[142-82-5]	17	%		
	n-Octane	[111-65-9]	17	%		
	n-Decane	[124-18-5]	12	%		
	n-Dodecane	[112-40-3]	5	%		
	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	12	%		
	Benzene	[71-43-2]	1	%		
	Toluene	[108-88-3]	9	%		
	Ethylbenzene	[100-41-4]	3	%		
	m-Xylene	[108-38-3]	3	%		
	o-Xylene	[95-47-6]	3	%		
	1,2,4-Trimethylbenzene	[95-63-6]	3	%		
	1,2,4,5-Tetramethylbenz ene	[95-93-2]	2	%		
	Naphthalene	[91-20-3]	1	%		