

## **1-Adamantyl Isocyanate**

Art. ID TRC-A220233-250MG  
Unit 250 mg  
Deliverydetails De Minimis /not restricted

### Description

Application Note: 1-Adamantyl Isocyanate is used to prepare reaction of nickel(III) diketiminato imido complex with CO.

| Text/Information | Analyte/Parameter      | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|-------------|---------------------|------|--------|--------|
|                  | 1-Adamantyl isocyanate | [4411-25-0] |                     |      |        |        |