

Adipic Acid-d4

Art. ID TRC-A291594-100MG
Unit 100 mg
Deliverydetails De Minimis

Description

Application Note: Hexanedioic-3,3,4,4-d4 Acid is the labeled version of A291590. Adipic Acid is primarily used in the synthesis of nylon. It has been used as a reagent in the solid-state polymerization of nylon analogs.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adipic Acid-d4	[121311-78-2]				