

Adipic Acid Propyl 2-Propyl Ester

Art. ID TRC-A291618-100MG
Unit 100 mg
Deliverydetails De Minimis

Description

Application Note: Adipic Acid Propyl 2-Propyl Ester is a derivative of Adipic Acid (A291590), which is primarily used in the synthesis of nylon. It has been used as a reagent in the solid-state polymerization of nylon analogs.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adipic Acid Propyl 2-Pr opyl Ester					