

L-Alanine-d3

Art. ID TRC-A481503-100MG
Unit 100 mg
Deliverydetails De Minimis /not restricted

Description

Application Note: L-Alanine-d3 is the labeled analogue of L-Alanine (A481500), which is used to make in-vivo measurement of glucose and alanine metabolism in studies of patients with diabetes. L-Alanine is a non-essential amino acid for human development. L-Alanine is one of the 20 amino acids encoded by the genetic code.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Alanine-3,3,3-d3	[63546-27-0]				