

3,7-Diamino-5-phenyl-5,10-dihydrophenazine

Art. ID TRC-A510990-100MG
Unit 100 mg
Deliverydetails De Minimis

Description

Application Note: 3,7-Diamino-5-phenyl-5,10-dihydrophenazine is a reduced form of Phenosafranin (P318020), which is a bacterial stain used in microscopy. 3,7-Diamino-5-phenyl-5,10-dihydrophenazine itself is patented for use as part of electrorecording material.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,7-Diamino-5-phenyl-5, 10-dihydrophenazine	[54668-92-7]				