

## **L-Alanyl-L-glutamine**

Art. ID TRC-A511370-25G  
Unit 25 g  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

### Description

Application Note: L-Alanyl-L-glutamine is a stable synthetic glutamine dipeptide. This dipeptide is used as an alternative source of amino acids in cell culture medium.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Alanyl-L-Glutamine	[39537-23-0]				