

2-Allyl-2-phenyl-4-pentenenitrile

Art. ID	TRC-A549215-10G
Unit	10 g
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Application Note: 2-Allyl-2-phenyl-4-pentenenitrile is an intermediate in synthesizing alpha-Phenyloxiranepropanenitrile (P335835), used in the preparation of CC-1065 and its analogs.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Allyl-2-phenyl-4-pentenenitrile	[28049-67-4]				