

DIN EN ISO 10301 - Halogenated VOCs, 1 µg/mL in Methanol

Art. ID ACCU-DINENISO-10301
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methylene chloride (Dichloromethane)	[75-09-2]	1	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	1	µg/mL		
	Carbon tetrachloride (Tetrachloromethane)	[56-23-5]	1	µg/mL		
	1,1-Dichloroethane	[75-34-3]	1	µg/mL		
	1,2-Dichloroethane	[107-06-2]	1	µg/mL		
	1,1,1-Trichloroethane	[71-55-6]	1	µg/mL		
	1,1,2-Trichloroethane	[79-00-5]	1	µg/mL		
	Trichloroethene	[79-01-6]	1	µg/mL		
	Tetrachloroethene	[127-18-4]	1	µg/mL		
	1,2-Dichloropropane	[78-87-5]	1	µg/mL		
	1,3-Dichloropropane	[142-28-9]	1	µg/mL		
	1,3-Dichloropropene (cis + trans)	[542-75-6]	1	µg/mL		
	Dibromomethane	[74-95-3]	1	µg/mL		
	Tribromoethene	[598-16-3]	1	µg/mL		
	Bromochloromethane	[74-97-5]	1	µg/mL		
	Bromodichloromethane	[75-27-4]	1	µg/mL		
	Dibromochloromethane	[124-48-1]	1	µg/mL		