

## **4-Aminophenethyl Alcohol**

Art. ID TRC-A585045-10G  
Unit 10 g  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

### Description

Application Note: 4-Aminophenethyl Alcohol, is a building block that can be used in the synthesis of various compounds, such as in synthesis of 4-aminostyrene, and in functionalization of graphene nanoplatelets.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Aminophenethyl alcoho	[104-10-9]				
	I					