

DL-2-Aminocaprylic Acid

Art. ID	TRC-A603185-5G
Unit	5 g
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: DL-2-Aminocaprylic Acid is used as a reagent in the preparation of 1,5-disubstituted-2-aminoimidazoles which has antibiotic activity against a variety of bacterial strains and even multi drug resistant isolates. DL-2-Aminocaprylic Acid is also used as a reactant in the solution phase parallel synthesis of fenbufen (F247200) and ethacrynic acid (E676000).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-2-Aminocaprylic acid	[644-90-6]				