

(trans-4-Aminocyclohexyl)acetic Acid Ethyl Ester

Art. ID	TRC-A604658-5G
Unit	5 g
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Application Note: (trans-4-Aminocyclohexyl)acetic Acid Ethyl Ester was used as a reagent in various reactions in the studies of modulator of the dopamine D2 receptor, D3 receptors for their uses in treatment of CNS disorders.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(trans-4-Aminocyclohexyl)acetic Acid Ethyl Ester	[76308-28-6]				