

5-Amino-6-nitroquinoline

Art. ID TRC-A618555-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: 5-Amino-6-nitroquinoline is a reagent used in the synthesis of benzodiazepine sulfonamide-based bombesin receptor 3 agonists. It is also used in the synthesis of quinoxalinediones.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------|--------------|---------------------|------|--------|--------|
| | 5-Amino-6-nitroquinolin e | [35975-00-9] | | | | |