

L-Arabinitol

Art. ID TRC-A735610-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: L-Arabinitol is the stereoisomer of D-Arabinitol (A735600), both used as diagnostic tools for the detection of fungal infections in neutropenic patients.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L(-)-Arabit	[7643-75-6]				