

Acid Red 114 (Technical Grade)

Art. ID	TRC-A795263-2.5G
Unit	2.5 g
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: Acid Red 114 is an azo dye used in diagnostics and textiles industry. It is also present in wastewater and textile effluents contributing to environmental hazards. Environmental toxin on US EPA Toxic Release Inventory list (TRI) list.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acid Red 114 (Technical Grade)	[6459-94-5]				