

4-Azidoaniline Hydrochloride

Art. ID TRC-A822020-2.5G
Unit 2.5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2926 (6.1) II

Description

Application Note: 4-Azidoaniline Hydrochloride is used to modify gelatin and bovine serum albumin in their photochemical immobilization on Ti-O films.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Azidoaniline hydrochl oride	[91159-79-4]				