

Bismuth ICP Standard, 1000 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-06N-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Bismuth (Bi) | [7440-69-9] | 1000 | µg/mL | | |