

## **Sulphide Ore Mill Tailings Sample**

Art. ID                      CANMET-RTS-3a  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

### **Description**

CANMET-RTS-3 was prepared from samples drawn by Noranda Inc. from the oxidized (vadose) zone and the saturated (unoxidized) zone of the Waite-Amulet tailings, located approximately 20 km north of Noranda, Quebec.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | S total           |             | 9,98 ± 0,26         | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 15,99 ± 0,46        | %    |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 2,20 ± 0,09         | %    |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 4,79 ± 0,17         | %    |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 2,45 ± 0,06         | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 2820 ± 90           | µg/g |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 1850 ± 80           | µg/g |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 146 ± 20            | µg/g |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 260 ± 16            | µg/g |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 9,1 ± 2,6           | µg/g |        |        |