

Potassium ICP Standard, 1000 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-43N-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Potassium (K) | [7440-09-7] | 1000 | µg/mL | | |