

**Tin ICP Standard, 10000 µg/mL in 2-5% Nitric Acid, tr Hydrofluoric acid**

Art. ID ACCU-ICP-63N-10X-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | Tin (Sn)          | [7440-31-5] | 10000               | µg/mL |        |        |