

**Copper ICP-MS Standard, 100 µg/mL in 2-5% Nitric Acid**

Art. ID                      ACCU-ICP-MS-15N-0.01X-1  
Unit                         100 mL  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | Copper (Cu)       | [7440-50-8] | 100                 | µg/mL |        |        |