

ICP-MS Calibration Standard 2 for Method 200.8, 10 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-MS-200.8-CAL2-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	µg/mL		
	Silver (Ag)	[7440-22-4]	10	µg/mL		